


## CHEMISTRY



## Inventory

S.No	Materials	Quantity
<b>Glassware</b>		
1	Beaker (Glass) – 250 ml	1
2	Dropper (1 ml)	2
3	Gas jar (plastic)	1
4	Test tubes	6
5	Glass rod	1
<b>Equipment</b>		
6	Aluminium wire	1 meter
7	Copper wire	1 meter
8	Deflagrating spoon	1
9	Electronic configuration model with marbles	1
10	Hammer	1
11	Spatula	2
12	Spirit lamp	1
13	Test tube stand	2
14	Burette stand	1
15	Plastic wire	1 meter
16	Wooden stick	1 meter
17	Blue litmus paper	1 packet
18	Red litmus paper	1 packet
19	1.5 V bulbs	5
20	1.5 V bulb holders	5
<b>Chemicals and Reagents</b>		
21	Copper sulphate	100 gm
22	Copper turnings	100 gm
23	Hydrochloric acid	500 ml
24	Iron fillings	100 gm
25	Iron sulphate	100 gm
26	Zinc metal	100 gm
27	Potassium permanganate	100 gm
28	Sulphur powder	50 gm
<b>Consumables</b>		
29	1.5 V Cell	2
30	A4 Sheets	1 bundle
31	Fevicol (125 gm)	1
32	Iron nails	10
33	Match box	1
34	Pencil	1
35	Thermocol balls	3 packets
36	Scale (long – 30 cm)	1

## Contents

S.NO	TOPIC	EXPERIMENT	PAGE NO
1	Electronic configuration	Electronic configuration	4
2	States of matter	Molecular arrangement in solid, liquid, gases state.	6
2.1	Dilution method		8
3	Metals and non metals	Malleability property of metals	10
3.1		Conductivity in metals	11
3.2		Reaction of metals	12
3.3		Reaction of non-metals	14
3.4		Reaction of acid with metals	15

## 1. Electronic configuration

### Aim

To explain the electronic configuration of elements according to Bohr's Model

### Materials Required

Electronic configuration board (model), marbles

### Procedure

#### Step 1

Take the electronic configuration board and count the holes of each different colour in each circle.

#### Step 2

Each circle represents a shell and assume that a hole in each shell can accommodate an electron.

#### Step 3

For example, consider sodium metal which has 11 electrons. Take 11 marbles and try to arrange the marbles in the electron configuration board. Begin by filling marbles in the innermost orbit as shown. Similarly we can try with different elements present in the periodic table.



**Observation**

In the periodic table, all elements have different number of electrons, and they are arranged in their respective shells.

**Inference**

- According to Bohr's model, there are shells around the nucleus, namely K, L, M, N ....
- The inner most shell is K and L, M, N are the next outer shells respectively.
- K-shell can hold 2 electrons.
- L-shell can hold 8 electrons
- M-shell can hold 18 electrons.
- N-shell can hold 32
- Electrons are filled in the respective shells accordingly.
- For Sodium, there are 2 electrons in K shell, 8 electrons in L shell and 1 electron in M shell
- Today, further studies have been made on the structure of atom. Bohr's model is only an approximate model.

## 2. States of matter

### Aim

To show molecular arrangement in solid, liquid and gas.

### Materials Required

A4 sheets, scale, pencil, thermocol balls and fevicol

### Procedure

#### Step1

Take A4 sheet and draw three equal squares (These three squares represent solid, liquid, and gaseous states respectively)

Step2 In the first square, stick the thermocol balls with the help of fevicol very closely. There should not be any gap between the balls. The first square represents the molecules in solid state.

Step3 Stick the thermocol balls in second square with small gaps between them. The second square represents the molecules in liquid state (due to random orientation of molecules).

Step4 In third square, stick the thermocol balls leaving big gaps between them as showed in the diagram. This square represents the molecules in gaseous state (due to more randomly arrangement of molecules).




### Observation

In the solid state the space between the molecules is very very less, in liquids the space between the molecules is more as compared to solids but in gases the space between the molecules is even more.

### Inference

In a solid state, the molecules are very closely packed because of the higher van der Waal's force of attraction between molecules; they cannot move freely and can only vibrate.

But in the case of liquids, the inter-molecular force of attraction is lower, so the molecules can move from one place to another within the volume of liquid.



In gaseous molecules, molecular arrangement is entirely different because the intermolecular force of attraction is very very weak, so the molecules are free to move from one place to another. They can spread out if not contained in space.

## 2.1 Dilution method

### Aim

To demonstrate the dilution of a solution.

### Materials required

Test tube, dropper, beaker, glass rod, Potassium permanganate, test tube stand

### Procedure

#### Step 1

Take five test tubes and label them as 1, 2, 3, 4 and 5.

#### Step 2

Take 10 ml of water in the first test tube and now add 3-4 potassium permanganate crystals into that test tube and stir it using a glass rod.

#### Step 3


Take 1 ml of solution from the first test tube and pour this solution into the second test tube which contains 9 ml of water.

#### Step 4

Take 1 ml of solution from the second test tube and pour this solution into the third test tube with 9 ml of water.

#### Step 5

Repeat the same procedure serially up to the fifth test tube as shown in the figure. Compare the intensity of colour among the test tubes.



### Observation

It is observed that the water in the last test tube is still colored but the intensity of colour is less when compared with other test tubes.

### Inference

When the potassium permanganate is dissolved in the first test tube, the concentration of  $\text{KMnO}_4$  is higher, so the higher colour of the solution. On



dilution, the concentration of  $\text{KMnO}_4$  decreases so intensity of the colour also decreases.

A solution is made by adding few crystals (x grams) of  $\text{KMnO}_4$  into 10 ml of water. Say the concentration is C.

In step 3, to get dilution, 1 ml of solution is added to 9 ml of water, concentration becomes 0.1 C

In step 4, it is diluted to .01C and finally to .001C in step 5. Clearly the intensity of colour decreases as solution becomes more and more dilute.

It is also understood from the activity that the minute particles of  $\text{KMnO}_4$  are separately dispersed in the solvent.

MAR 201

### 3. Malleability property of metals

**Aim**

To demonstrate the malleability property of metals and non-malleability or brittleness of non metals

**Materials Required**

Aluminium wire, iron nail, copper wire, magnesium wire, charcoal, graphite and hammer.

**Procedure****Step 1**

Take iron nail, keep it on the floor and hit the iron nail using hammer.

**Step 2**

Similarly repeat for the rest of the materials.

What do you observe?

**Observation**

Metal objects get flattened upon hammering while non-metal breaks.

**Inference**

Metals are malleable (Ex: aluminum wire, iron nail, copper wire etc) where as non-metals (Ex: graphite and charcoal) are brittle.

### 3.1 Conductivity in metals

#### Aim

To demonstrate the conductivity of metals and graphite.


#### Materials Required

Aluminium wire, iron nail, copper wire, magnesium wire, pencil, piece of glass, wood, dry cell, a bulb (low voltage) with base and terminals. (bulb holders).

#### Procedure

##### Step 1

Connect the cell, bulb and copper connecting wires as shown in the figure. Place the aluminum wire in contact, with free ends of the connecting wires. The bulb glows.



##### Step 2

Similarly repeat the experiment with other materials provided.

#### Observation

The bulb glows when metals and graphite are connected to the cell, whereas for non metals like wood, plastic and glass, the bulb doesn't glow.

#### Inference

Metals and graphite are good conductors of electricity due to presence of free electrons in their outer most orbits. In case of wood, plastic and glass, there are no free electrons, hence they are insulators in nature.



### 3.2 Reaction of metals

#### Aim

To prove that iron is more reactive than copper

#### Materials Required

Iron fillings, copper turnings, copper sulphate solution, iron sulphate solution, test tubes, test tube stand and dropper.



#### Procedure

##### Step 1

Take few gm of iron fillings in a test tube and add half test tube full of copper sulphate solution with the help of a dropper (test tube –A).

##### Step 2

Take few gm of copper turnings in a test tube and add half test tube full of iron sulphate solution with the help of another dropper (test tube B).

### Observation

In test tube A, iron fillings turns to brown and the solution turns to pale green. No reaction is observed in the test tube B


### Inference

In test tube A, iron fillings turn to brown due to the deposition of copper. Because  $\text{Cu}^{2+}$  ions are easily reduced to copper atom by accepting electrons from iron and the solution turns to pale green and the formation of iron sulphate.


But in case of test tube B, no reaction takes place with respect to copper fillings and iron sulphate solution because copper in reduced state, will not oxidise easily. Iron therefore is more reactive than copper as shown in the reaction.

The reaction can be written as

#### Test tube A



i.e.



#### Test tube B

No reaction

### 3.3 Reaction of non-metals

#### Aim

To prepare a non metallic oxide and to test its solution using litmus paper

#### Materials Required

Sulphur powder, deflagrating spoon, spatula, gas jar, water, spirit lamp, match box and blue litmus paper

#### Procedure

##### Step 1

Place a little sulphur powder in the deflagrating spoon and heat it. When it starts to burn, place it in the gas jar.

##### Step 2

Pour about 50 ml of water in to the gas jar.

##### Step 3



Dip the blue litmus paper in the water. Observe the colour of litmus paper.

#### Observation

Blue litmus paper turns red.

#### Inference

Sulphur burns in presence oxygen and forms sulphur dioxide. When water is added to sulphur di oxide, it forms sulphurous acid.



Note: Avoid sulphur contact with eyes, mouth or mucous membrane because it is corrosive



### 3.4 Reaction between acid and metal

#### Aim

To demonstrate the displacement reaction between acid and metal

#### Materials Required

Test tubes, burette stand, zinc metal pieces, dilute hydrochloric acid and match box

#### Procedure

##### Step 1

Take a test tube and clamp it to the burette stand. Put few zinc metal pieces in to the test tube.

##### Step 2

By means of dropper add dilute hydrochloric acid into the test tube

##### Step 3

Light a matchstick and bring it near the mouth of the test tube.

#### Observation

When the match stick is placed at the mouth of the test tube, the evolved gas catches a fire and gives a bursting sound while the match stick is extinguished.

#### Inference

Zinc metal is highly reactive with hydrochloric acid, it displaces hydrogen to form zinc chloride with the release of hydrogen gas.

